

nivan

NIVAN CARE

# PRODUCT CATALOG



PRODUCT CATALOG



WWW.NIVANCARE.COM

# TABLE OF CONTENT

## Tecar(RF) Cream

03

Specially formulated for use with Tecar(RF) therapy devices, Tecar(RF) Cream ensures optimal conductivity, smooth application, and enhanced comfort during treatment. Its balanced texture allows easy glide while supporting effective energy transfer to the targeted area.

## Ultrasonic Gel

06

Ultrasonic Gel is designed for diagnostic and therapeutic ultrasound procedures, providing excellent transmission of ultrasonic waves. Its non-irritating, water-based formula ensures clarity, comfort, and consistent performance during each session.



# TECAR(RF) CREAM



Specialized Tecar (RF) Creams are formulated for a variety of therapeutic and aesthetic purposes:

**Normal Cream:** The best choice for frequent and general use, offering a smooth texture suitable for all skin types.

**Body Plus Cream:** Designed for body contouring and slimming treatments, enriched with anti-cellulite and fat-burning ingredients to enhance skin tone and texture.

**Face & Body Cream:** Ideal for slimming and beauty treatments in the face and neck area, containing lifting, rejuvenating, moisturizing, and fat-burning agents. Its smoothing components help reduce skin redness during treatment.





# TECAR(RF) CREAM



Tecar Creams are available in two packaging options:

Jar – 1000 ml, ideal for very frequent and high-volume clinical use.

Tube – 250 ml, more hygienic and easier to use.



nivan

# NEW PRODUCT

Tecar PT Plus Cream, specially developed for physiotherapy applications.

With its advanced formula, it effectively supports the treatment process and provides extra assistance in reducing pain, ensuring a more comfortable therapy experience.





# ULTRASONIC GEL



Ultrasonic Gel is specially designed for diagnostic and therapeutic use with ultrasound devices.

This high-quality, thick gel contains glycerin, ensuring optimal transmission of ultrasonic waves. Its water-based, smooth, and non-irritating formula enhances patient comfort while improving clarity and precision in imaging or treatment.

# ULTRASONIC GEL

Ultrasonic Gels are available in two packaging options:

5-liter gallon – ideal for high-volume and frequent clinical use.

500 ml bottle – more hygienic and easier to use.







## Economic Ultrasonic Gels

Our new product, Economic Ultrasonic Gels, is designed for wide applications, offering good quality at an affordable price. They are available in three packaging options:

5-liter gallon – ideal for very high-volume clinical use.

1-liter bottle – a balanced choice for medium-level consumption.

260 ml bottle – more hygienic and easier to use, perfect for personal or limited applications.

# NEW PRODUCT



**5L**

**epulse<sub>gel</sub>**  
Ultrasound

**Diagnostic Sonography, Shockwave Therapy  
Ultrasound Therapy, EEG, EMG, ECG, Laser Therapy**

• No color • Water soluble • Non-irritating • Dermatologically tested  
 • External use only • No formaldehyde • Electrically conductive  
 • No fragrance • Hypoallergenic • Non-greasy

سازمان غذا و دارو  
جمهوری اسلامی ایران  
۳۱۴۱۶۲۱۱

3 Years after PROD.DATE

Made in Iran  
تولید شده در ایران

**NOVINMED**  
A Lifelong Friend

No. 169, Electrical & Electronics Industrial Zone of  
 Khomeinshahr, Isfahan, Islamic Republic of Iran  
[www.novinmed.com](http://www.novinmed.com) [info@novinmed.com](mailto:info@novinmed.com)  
 Tel: (+98)2141223

nivan





# Tecar (RF) cream

|                | Usage Area  | Advantage                                                         | Pakaging                   |
|----------------|-------------|-------------------------------------------------------------------|----------------------------|
| NORMAL         | Body        | public use<br>Economic                                            | Jar 1000 ml<br>Tube 250 ml |
| BODY+          | Body        | Anti-cellulit<br>Brightening<br>Slimming                          | Jar 1000 ml<br>Tube 250 ml |
| FACE &<br>BODY | Face & Body | Brightening<br>Slimming<br>Anti-wrinkle<br>Hydrating<br>Smoothing | Jar 1000 ml<br>Tube 250 ml |
| PT+            | Body        | Relief                                                            | Jar 1000 ml<br>Tube 250 ml |

# Ultrasonic Gel

|          | Usage Area  | Advantage                                          | Pakaging                                      |
|----------|-------------|----------------------------------------------------|-----------------------------------------------|
| EPULSE + | Face & Body | concentrated gel<br>Hypoallergenic<br>Long lasting | Bottle 500 ml<br>Gallon 5 L                   |
| EPULSE   | Body        | Public use<br>Economic                             | Bottle 260 ml<br>Bottle 1000 ml<br>Gallon 5 L |

# CONTACT US



+989367298301  
+989361848301



nivancare



[www.nivancare.com](http://www.nivancare.com)



[Nivan@Nivancare.com](mailto:Nivan@Nivancare.com)



Unit 7, Golestan Building  
Next to Park Complex  
Chaharbagh Bala St.  
Isfahan  
Iran

nivan

